

Technical Publication

DR-1004R13

LV-DRAC

Low Voltage - Digital Rotating Anode Controller

HF Series Generators

REVISION HISTORY

REVISION	DATE	REASON FOR CHANGE
6	MAR 1, 2007	New Tubes in Table-2 (soft V10R2.04)
7	MAY 16, 2007	New Tubes in Table-2 (soft V10R3.1)
8	OCT 16, 2007	New Tubes in Table-2 (soft V10R3.3)
9	JAN 22, 2008	New Tubes in Table-2 (soft V10R3.4)
10	MAY 12, 2008	New Tubes in Table-2 (soft V10R3.5)
11	JUL 01, 2009	New Tubes and schematics (soft V10R3.6)
12	MAR 17, 2011	New Tubes (soft V10R3.7)
13	MAY 14, 2012	New Tubes (soft V10R3.8.1)

This Document is the English original version, edited and supplied by the manufacturer.

The Revision state of this Document is indicated in the code number shown at the bottom of this page.

ADVISORY SYMBOLS

The following advisory symbols will be used throughout this manual. Their application and meaning are described below.



DANGERS ADVISE OF CONDITIONS OR SITUATIONS THAT IF NOT HEDED OR AVOIDED WILL CAUSE SERIOUS PERSONAL INJURY OR DEATH.



ADVISE OF CONDITIONS OR SITUATIONS THAT IF NOT HEDED OR AVOIDED COULD CAUSE SERIOUS PERSONAL INJURY, OR CATASTROPHIC DAMAGE OF EQUIPMENT OR DATA.



Advise of conditions or situations that if not heeded or avoided could cause personal injury or damage to equipment or data.

Note 

Alert readers to pertinent facts and conditions. Notes represent information that is important to know but which do not necessarily relate to possible injury or damage to equipment.

TABLE OF CONTENTS

Section		Page
1	CONFIGURATION	1
	1.1 LV-DRAC Status Diagram	10
	1.2 Self-running Mode	10
2	SCHEMATICS	11

SECTION 1 CONFIGURATION

Note 

This document ONLY applies to High Speed Generators supplied with the LV-DRAC (Low Voltage - Digital Rotating Anode Controller).

The LV-DRAC is a High Speed Rotor Controller located on the Generator Cabinet. Its configuration is made through dip switches 3243SW1 to 3243SW4 on the Control DRAC Board.



BEFORE MANIPULATING THE LV-DRAC, MAKE SURE THAT THE INPUT LINE IS DISCONNECTED AND THE CAPACITOR BANK IS PROPERLY DISCHARGED. WAIT UNTIL LEDS DL7 AND DL8 ON THE CONTROL DRAC BOARD ARE OFF.

The Function of these switches is the following:

3243SW1-x	FUNCTION
1	MINIMUM TIME FOR READY (DELAY)
2	
3	
4	TUBE TYPE SELECTION FOR TUBE-1 (Also refer to switch position 3243SW3-5)
5	
6	
7	
8	

3243SW2-x	FUNCTION
1	FLUORO HOLD-OVER TIME
2	
3	
4	TUBE TYPE SELECTION FOR TUBE-2 (Also refer to switch position 3243SW4-5)
5	
6	
7	
8	

3243SW3-x	FUNCTION
1	LOW SPEED ROTATION BRAKE
2	NOT USED*
3	SPOT FILM HIGH SPEED START
4	NOT USED*
5	ON = WHEN THE TUBE CONFIGURED AS TUBE-1 IS LISTED ON TUBE TABLE-1
	OFF = WHEN THE TUBE CONFIGURED AS TUBE-1 IS LISTED ON TUBE TABLE-2
6	NOT USED*
7	
8	
* Note: Set switches not used in "ON" position.	

3243SW4-x	FUNCTION
1	SPOT FILM HOLD TIME
2	
3	
4	
5	ON = WHEN THE TUBE CONFIGURED AS TUBE-2 IS LISTED ON TUBE TABLE-1
	OFF = WHEN THE TUBE CONFIGURED AS TUBE-2 IS LISTED ON TUBE TABLE-2
6	ON = SELF RUNNING ENABLED (Refer to section 1.2)
	OFF = SELF RUNNING DISABLED (Refer to section 1.2)
7	PROTECTIONS - ERRORS
8	DC BRAKE
* Note: Set switches not used in "ON" position.	

Note 

Set dip switches in accordance with the following tables for Conventional or Mammographic use.

TUBE SELECTION - TUBE TABLE-1 for CONVENTIONAL USE

(Switch 3243SW3-5 = ON for Tube-1) (Switch 3243SW4-5 = ON for Tube-2)

Pos	TUBE FAMILY (Stator - ∅ Anode)	3243SW1-x (TUBE-1)					3243SW2-x (TUBE-2)				
		4	5	6	7	8	4	5	6	7	8
0	GE MAXIRAY 75	ON	ON	ON	ON	ON	ON	ON	ON	ON	ON
1	GE MAXIRAY 100	OFF	ON	ON	ON	ON	OFF	ON	ON	ON	ON
2	VARIAN G-1582	ON	OFF	ON	ON	ON	ON	OFF	ON	ON	ON
3	VARIAN RAD-XX or EUREKA RAD-XX (3" Anode, Emerald or Diamond Housing, R Stator, 20/50 ohm, 300 KHU)	OFF	OFF	ON	ON	ON	OFF	OFF	ON	ON	ON
	VARIAN A132 (3" Anode, 300 KHU)										
4	VARIAN RAD-XX or EUREKA RAD-XX (4" Anode, Sapphire Housing, R Stator (16/50 ohm or 20/50 ohm), 400 KHU)	ON	ON	OFF	ON	ON	ON	ON	OFF	ON	ON
	VARIAN RAD-56 (4" Anode, Sapphire Housing, R Stator (20/50 ohm), 400 KHU)										
	VARIAN A256 / A282 / A292 (4" Anode, B130 or B150 Housing, R Stator (16/50 ohm or 20/50 ohm), 400 KHU)										
	TOSHIBA 7254X / 7255X / 7823 (4" Anode, XS-RB Stator (20/38 ohm), 400 KHU)										
5	VARIAN GS-2075 with B220 Housing	OFF	ON	OFF	ON	ON	OFF	ON	OFF	ON	ON
6	SIEMENS 100 L BIANGULIX SIEMENS 100 L OPTILIX	ON	OFF	OFF	ON	ON	ON	OFF	OFF	ON	ON
7	PHILIPS SUPEROTALIX	OFF	OFF	OFF	ON	ON	OFF	OFF	OFF	ON	ON
8	DUNLEE PX-1302	ON	ON	ON	OFF	ON	ON	ON	ON	OFF	ON
9	CGR STATORIX 260 with RSN 742 / MSN 740	OFF	ON	ON	OFF	ON	OFF	ON	ON	OFF	ON
10	CGR STATORIX 550 with RSN 722 / MSN 722	ON	OFF	ON	OFF	ON	ON	OFF	ON	OFF	ON
11	SIEMENS 100 G	OFF	OFF	ON	OFF	ON	OFF	OFF	ON	OFF	ON
12	SIEMENS RH 150/100	ON	ON	OFF	OFF	ON	ON	ON	OFF	OFF	ON
13	COMET DO-10 (4" Anode)	OFF	ON	OFF	OFF	ON	OFF	ON	OFF	OFF	ON
14	COMET DO-9 (4" Anode)	ON	OFF	OFF	OFF	ON	ON	OFF	OFF	OFF	ON
15	CGR STATORIX 240 with MN 620 / M 641	OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF	ON
16	PICKER PX-307	ON	ON	ON	ON	OFF	ON	ON	ON	ON	OFF
17	HANGZHOU XD-52-30, 50/150	OFF	ON	ON	ON	OFF	OFF	ON	ON	ON	OFF
18	DUNLEE PX-1312	ON	OFF	ON	ON	OFF	ON	OFF	ON	ON	OFF
19	DUNLEE PX-1456 / PX-1436 / PX-1400 (Q Stator)	OFF	OFF	ON	ON	OFF	OFF	OFF	ON	ON	OFF
20	COMET DX-1000	ON	ON	OFF	ON	OFF	ON	ON	OFF	ON	OFF
21	COMET DX-81 / COMET XSTAR-14	OFF	ON	OFF	ON	OFF	OFF	ON	OFF	ON	OFF
22	PHILIPS ROTALIX 350/351	ON	OFF	OFF	ON	OFF	ON	OFF	OFF	ON	OFF
23		OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF	ON	OFF
24		ON	ON	ON	OFF	OFF	ON	ON	ON	OFF	OFF
25	VARIAN with B160 / B165 Housing	OFF	ON	ON	OFF	OFF	OFF	ON	ON	OFF	OFF
26	TOSHIBA E7252X (XS-AL Stator)	ON	OFF	ON	OFF	OFF	ON	OFF	ON	OFF	OFF
27	COMET DX 700 HS	OFF	OFF	ON	OFF	OFF	OFF	OFF	ON	OFF	OFF
28	VARIAN A182 / A192 / A196 / A197 (with B130 Housing)	ON	ON	OFF	OFF	OFF	ON	ON	OFF	OFF	OFF
29	VARIAN G 1592	OFF	ON	OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF
30	PICKER PX400P + PX457P	ON	OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF
31	JUGORENDGEN RX 150/30 - 50	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF

HF Series Generators

LV-DRAC

TUBE SELECTION - TUBE TABLE-2 for CONVENTIONAL USE

(Switch 3243SW3-5 = OFF for Tube-1) (Switch 3243SW4-5 = OFF for Tube-2)

Pos	TUBE FAMILY (Stator - ∅ Anode)	3243SW1-x (TUBE-1)					3243SW2-x (TUBE-2)				
		4	5	6	7	8	4	5	6	7	8
0	VARIAN GS-20711 (B220H Housing, R Stator)	ON	ON	ON	ON	ON	ON	ON	ON	ON	ON
1	SIEMENS 100 L-W	OFF	ON	ON	ON	ON	OFF	ON	ON	ON	ON
2	SIEMENS SV 125/15/82	ON	OFF	ON	ON	ON	ON	OFF	ON	ON	ON
3	IAE RTM90 / IAE RTM92 (C52 / C352 Housing)	OFF	OFF	ON	ON	ON	OFF	OFF	ON	ON	ON
4	IAE RTM101HS / IAE RTM102HS (C100 Housing)	ON	ON	OFF	ON	ON	ON	ON	OFF	ON	ON
5	DUNLEE DA 1036 - DU 404	OFF	ON	OFF	ON	ON	OFF	ON	OFF	ON	ON
6	DUNLEE DA 1094 - DU 694 / DU 692	ON	OFF	OFF	ON	ON	ON	OFF	OFF	ON	ON
7	DUNLEE DR 1817 / 1825 (High Impedance Stator)	OFF	OFF	OFF	ON	ON	OFF	OFF	OFF	ON	ON
8	CGR STATORIX 240 / 260 with MSN 740 / MSN 742	ON	ON	ON	OFF	ON	ON	ON	ON	OFF	ON
9	IAE RTC 700 HS (C52S or C100XT Housing)	OFF	ON	ON	OFF	ON	OFF	ON	ON	OFF	ON
10	VARIAN GS 30711	ON	OFF	ON	OFF	ON	ON	OFF	ON	OFF	ON
11	VARIAN G292 / G294 / G256 / G297 (B130 Housing)	OFF	OFF	ON	OFF	ON	OFF	OFF	ON	OFF	ON
12	IAE RTM101HS / RTM102HS (C52S Housing)	ON	ON	OFF	OFF	ON	ON	ON	OFF	OFF	ON
13	IAE RTC1000HS (C52S Housing)	OFF	ON	OFF	OFF	ON	OFF	ON	OFF	OFF	ON
14	TOSHIBA ROTANODE E7100X (4" Anode, 9/28 ohm, RTM, 300 KHU)	ON	OFF	OFF	OFF	ON	ON	OFF	OFF	OFF	ON
15	DUNLEE DA 1083 / PX 1483 (DU 404 / DA10 Housing, C Stator)	OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF	ON
16	VARIAN A277 / A278 (B130H and B150H Housing, R Stator)	ON	ON	ON	ON	OFF	ON	ON	ON	ON	OFF
17	VARIAN G1092 (B160 / B165 Housing, R Stator)	OFF	ON	ON	ON	OFF	OFF	ON	ON	ON	OFF
18	DUNLEE PX1473Q - DU 604	ON	OFF	ON	ON	OFF	ON	OFF	ON	ON	OFF
19	TOSHIBA E7260DX	OFF	OFF	ON	ON	OFF	OFF	OFF	ON	ON	OFF
20	VARIAN HE100 (B180 / B185 Housing, R Stator)	ON	ON	OFF	ON	OFF	ON	ON	OFF	ON	OFF
21	IAE RTC 137 (CT 180 Housing)	OFF	ON	OFF	ON	OFF	OFF	ON	OFF	ON	OFF
22	VARIAN A102 (B100 Housing)	ON	OFF	OFF	ON	OFF	ON	OFF	OFF	ON	OFF
23	IAE X40 (C352 Housing (only Low Speed)) IAE RTM78HS (C352 Housing)	OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF	ON	OFF
24	VARIAN RAD-92 / RAD-94 (4" Anode, Sapphire Housing, R Stator, 20/50 ohm, 600 KHU) VARIAN SG296B (B199-R Housing) IAE RTC600 HS (4" Anode, C100 Housing, 20/40 ohm, 600 KHU)	ON	ON	ON	OFF	OFF	ON	ON	ON	OFF	OFF
25	TOSHIBA ROTANODE E7867X / E7869X	OFF	ON	ON	OFF	OFF	OFF	ON	ON	OFF	OFF
26	TOSHIBA ROTANODE E7239X / E7240X / E7242X / E7865X	ON	OFF	ON	OFF	OFF	ON	OFF	ON	OFF	OFF
27	TOSHIBA ROTANODE E7864X (4" Anode, 9/28 ohm, 400 KHU)	OFF	OFF	ON	OFF	OFF	OFF	OFF	ON	OFF	OFF
28	SIEMENS RAY-14_1 (SV150/33/78S Insert)	ON	ON	OFF	OFF	OFF	ON	ON	OFF	OFF	OFF
29	TOSHIBA ROTANODE E7252X (XS-R Stator) TOSHIBA ROTANODE E7886X (27/58 ohm Stator)	OFF	ON	OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF
30	ATX 10WX-692 / 1463 (AT1692) (4" Anode, C Stator, 47/60 ohm, 600 KHU)	ON	OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF
31	PHILIPS SRM 1551 ROT500	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF

TUBE SELECTION - TUBE TABLE-1 for MAMMOGRAPHIC USE

(Switch 3243SW3-5 = ON for Tube-1) (Switch 3243SW4-5 = ON for Tube-2)

Pos	TUBE FAMILY (Stator - ∅ Anode)	3243SW1-x (TUBE-1)					3243SW2-x (TUBE-2)				
		4	5	6	7	8	4	5	6	7	8
0		ON	ON	ON	ON	ON	ON	ON	ON	ON	ON
1		OFF	ON	ON	ON	ON	OFF	ON	ON	ON	ON
2		ON	OFF	ON	ON	ON	ON	OFF	ON	ON	ON
3	IAE XM1016T (3" Anode, 20/40 ohm, 300 KHU) (MAMMO)	OFF	OFF	ON	ON	ON	OFF	OFF	ON	ON	ON
4		ON	ON	OFF	ON	ON	ON	ON	OFF	ON	ON
5	CGR STATORIX M50 (MAMMO)	OFF	ON	OFF	ON	ON	OFF	ON	OFF	ON	ON
6		ON	OFF	OFF	ON	ON	ON	OFF	OFF	ON	ON
7		OFF	OFF	OFF	ON	ON	OFF	OFF	OFF	ON	ON
8		ON	ON	ON	OFF	ON	ON	ON	ON	OFF	ON
9		OFF	ON	ON	OFF	ON	OFF	ON	ON	OFF	ON
10		ON	OFF	ON	OFF	ON	ON	OFF	ON	OFF	ON
11		OFF	OFF	ON	OFF	ON	OFF	OFF	ON	OFF	ON
12		ON	ON	OFF	OFF	ON	ON	ON	OFF	OFF	ON
13		OFF	ON	OFF	OFF	ON	OFF	ON	OFF	OFF	ON
14		ON	OFF	OFF	OFF	ON	ON	OFF	OFF	OFF	ON
15		OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF	ON
16		ON	ON	ON	ON	OFF	ON	ON	ON	ON	OFF
17		OFF	ON	ON	ON	OFF	OFF	ON	ON	ON	OFF
18		ON	OFF	ON	ON	OFF	ON	OFF	ON	ON	OFF
19		OFF	OFF	ON	ON	OFF	OFF	OFF	ON	ON	OFF
20		ON	ON	OFF	ON	OFF	ON	ON	OFF	ON	OFF
21		OFF	ON	OFF	ON	OFF	OFF	ON	OFF	ON	OFF
22		ON	OFF	OFF	ON	OFF	ON	OFF	OFF	ON	OFF
23	COMET MOS-50 (MAMMO)	OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF	ON	OFF
24	VARIAN M113 with B115 Housing, R Stator (MAMMO) VARIAN IM113T (MAMMO) VARIAN VRAD 85S with B112 Housing (3", 16/50 ohm, 300 KHU) (MAMMO)	ON	ON	ON	OFF	OFF	ON	ON	ON	OFF	OFF
25		OFF	ON	ON	OFF	OFF	OFF	ON	ON	OFF	OFF
26		ON	OFF	ON	OFF	OFF	ON	OFF	ON	OFF	OFF
27		OFF	OFF	ON	OFF	OFF	OFF	OFF	ON	OFF	OFF
28		ON	ON	OFF	OFF	OFF	ON	ON	OFF	OFF	OFF
29		OFF	ON	OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF
30		ON	OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF
31		OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF

HF Series Generators

LV-DRAC

TUBE SELECTION - TUBE TABLE-2 for MAMMOGRAPHIC USE

(Switch 3243SW3-5 = OFF for Tube-1) (Switch 3243SW4-5 = OFF for Tube-2)

Pos	TUBE FAMILY (Stator - ∅ Anode)	3243SW1-x (TUBE-1)					3243SW2-x (TUBE-2)				
		4	5	6	7	8	4	5	6	7	8
0		ON	ON	ON	ON	ON	ON	ON	ON	ON	ON
1		OFF	ON	ON	ON	ON	OFF	ON	ON	ON	ON
2		ON	OFF	ON	ON	ON	ON	OFF	ON	ON	ON
3		OFF	OFF	ON	ON	ON	OFF	OFF	ON	ON	ON
4		ON	ON	OFF	ON	ON	ON	ON	OFF	ON	ON
5		OFF	ON	OFF	ON	ON	OFF	ON	OFF	ON	ON
6		ON	OFF	OFF	ON	ON	ON	OFF	OFF	ON	ON
7		OFF	OFF	OFF	ON	ON	OFF	OFF	OFF	ON	ON
8		ON	ON	ON	OFF	ON	ON	ON	ON	OFF	ON
9		OFF	ON	ON	OFF	ON	OFF	ON	ON	OFF	ON
10		ON	OFF	ON	OFF	ON	ON	OFF	ON	OFF	ON
11		OFF	OFF	ON	OFF	ON	OFF	OFF	ON	OFF	ON
12		ON	ON	OFF	OFF	ON	ON	ON	OFF	OFF	ON
13		OFF	ON	OFF	OFF	ON	OFF	ON	OFF	OFF	ON
14		ON	OFF	OFF	OFF	ON	ON	OFF	OFF	OFF	ON
15		OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF	ON
16		ON	ON	ON	ON	OFF	ON	ON	ON	ON	OFF
17		OFF	ON	ON	ON	OFF	OFF	ON	ON	ON	OFF
18		ON	OFF	ON	ON	OFF	ON	OFF	ON	ON	OFF
19		OFF	OFF	ON	ON	OFF	OFF	OFF	ON	ON	OFF
20		ON	ON	OFF	ON	OFF	ON	ON	OFF	ON	OFF
21		OFF	ON	OFF	ON	OFF	OFF	ON	OFF	ON	OFF
22		ON	OFF	OFF	ON	OFF	ON	OFF	OFF	ON	OFF
23		OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF	ON	OFF
24	VARIAN RAD 70 (MAM-RAD 105 H Housing, R Stator) (MAMMO) VARIAN RAD 70B (MAM-RAD 105 H Housing, R Stator) (MAMMO) VARIAN RAD 70SP (MAM-RAD 100 H Housing, R Stator) (MAMMO)	ON	ON	ON	OFF	OFF	ON	ON	ON	OFF	OFF
25		OFF	ON	ON	OFF	OFF	OFF	ON	ON	OFF	OFF
26		ON	OFF	ON	OFF	OFF	ON	OFF	ON	OFF	OFF
27		OFF	OFF	ON	OFF	OFF	OFF	OFF	ON	OFF	OFF
28		ON	ON	OFF	OFF	OFF	ON	ON	OFF	OFF	OFF
29		OFF	ON	OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF
30		ON	OFF	OFF	OFF	OFF	ON	OFF	OFF	OFF	OFF
31		OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF	OFF

MINIMUM TIME FOR READY (DELAY)

MINIMUM TIME (SECONDS)	3243SW1-x		
	1	2	3
0 *	ON	ON	ON
0.6	OFF	ON	ON
0.8	ON	OFF	ON
1.0	OFF	OFF	ON
1.2	ON	ON	OFF
1.5	OFF	ON	OFF
2.0	ON	OFF	OFF
3.0	OFF	OFF	OFF

** Note: It is recommended to set these switches at "0 seconds".*

FLUORO HOLD TIME

FLUORO HOLD TIME (SECONDS)	3243SW2-x		
	1	2	3
0	ON	ON	ON
10	OFF	ON	ON
20	ON	OFF	ON
30	OFF	OFF	ON
40	ON	ON	OFF
60	OFF	ON	OFF
90	ON	OFF	OFF
120	OFF	OFF	OFF

SPOT FILM HOLD TIME

SPOT FILM HOLD TIME		3243SW4-x			
		1	2	3	4
SECONDS	0	ON	ON	ON	ON
	5	OFF	ON	ON	ON
	10	ON	OFF	ON	ON
	15	OFF	OFF	ON	ON
	20	ON	ON	OFF	ON
	30	OFF	ON	OFF	ON
	45	ON	OFF	OFF	ON
	60	OFF	OFF	OFF	ON
MINUTES	2	ON	ON	ON	OFF
	5	OFF	ON	ON	OFF
	10	ON	OFF	ON	OFF
	15	OFF	OFF	ON	OFF
	20	ON	ON	OFF	OFF
	25	OFF	ON	OFF	OFF
	30	ON	OFF	OFF	OFF
	40	OFF	OFF	OFF	OFF

SPOT FILM HIGH SPEED START

SPOT FILM START MODE	3243SW3-3
SPOT FILM ALWAYS STARTS IN HIGH SPEED (10000 RPM)	ON
SPOT FILM STARTS AT REQUIRED SPEED	OFF

LOW SPEED ROTATION BRAKE

LOW SPEED ROTATION BRAKE MODE (3300 RPM TO 0 RPM)	3243SW3-1
INHIBITED	ON
ACTIVATED	OFF

DC BRAKE (FOR LOW AND HIGH SPEED)

DC BRAKE	3243SW4-8
ACTIVATED	ON
INHIBITED	OFF


PROTECTIONS - ERRORS

PROTECTIONS - ERRORS	3243SW4-7
INHIBITED (LED DL1 on the Control DRAC Board is always illuminated)	ON
ACTIVATED (Normal Operation mode) (LED DL1 on the Control DRAC Board is non-illuminated indicating that the DRAC is working properly).	OFF
<i>Note: The Error Codes are shown on the Console Display preceded by the letter "E" (a.e. E51) (For Error Codes refer to the Troubleshooting document).</i>	



WHEN SWITCH 3243SW4-7 IS IN THE "ON" POSITION, THE ERROR PROTECTIONS OF THE LV-DRAC ARE INHIBITED.

1.1 LV-DRAC STATUS DIAGRAM




1.2 SELF-RUNNING MODE

The continuous starting and braking of the anode produces an overheating in the X-ray Tube Stator (a.e. during calibration procedure). The “Self-running mode” avoids this overheating.

In the “Self-running mode” (switch 3243SW4-6 in the “ON” position), the anode remains running for approximately 1 minute when rotates at 3300 rpm or 10000 rpm.

The LV-DRAC enters in the “Self-running mode” when “Preparation” is selected three consecutive times **only from the Console buttons or Handswitch**, and the time between two of the consecutive accelerations is shorter than 30 seconds (refer to illustration below).




tsr (time for self-running mode) = approx. 1 minute at 3300 rpm or 10000 rpm.


SECTION 2 SCHEMATICS

Sch. No.	Scheme
	LV-DRAC MODULE BLOCK DIAGRAM
A3243-04	CONTROL DRAC PCB
A3240-05	INTERFACE DRAC PCB
A3109-01	CLAMPING PCB

This page intentionally left blank.



LV-DRAC Module Block Diagram



				NAME	DATE	SHEET / OF	A3243-04			
				DRAWING	F. Diaz	12/06/08				
				REVISED	A. Diaz	12/06/08				
F	NC 11/185	F. Diaz	30/05/11							
E	NC 09/204	F. Diaz	14/07/09							
D	NC 09/151	F. Diaz	03/05/09							
REV	DESCRIPTION	ISSUED BY	DATE							




CONTROL DRAC


PTR1

1/3

F E D ← REV




REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3243-04					
F	NC 11/185	F. Diaz	30/05/11	DRAWING	F. Diaz	12/06/08	F E D ← REV					
E	NC 09/204	F. Diaz	14/07/09	REVISED	A. Diaz	12/06/08	2/3					
D	NC 09/151	F. Diaz	03/05/09								CONTROL DRAC	




NOTE: 4.5 V < > 37 Amp. pick for I. PRINC or I. AUX

T1 IF: T1 to T2 Commutation
 T2 IF: Not Used
 T3 IF: Break Control

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3243-04				
				DRAWING	F. Diaz	12/06/08					
				REVISED	A. Diaz	12/06/08					
F	NC 11/185	F. Diaz	30/05/11								
E	NC 09/204	F. Diaz	14/07/09								
D	NC 09/151	F. Diaz	03/05/09								


CONTROL DRAC



PCB: 3240-05D

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3240-05						
G	NC 09/301	F. Diaz	16/11/09				A3240-05						
F	NC 09/178	F. Diaz	09/06/09	DRAWING	F. GARCIA	27-09-97	A3240-05						
E	CN 01/075	F. GARCIA	03-05-01	REVISED	A. DIAZ	27-09-97	A3240-05						
C, D	CN 98/090, 99/94	F. GARCIA	15-02-00				A3240-05						
H	NC 10/252	F. Diaz	28/05/10				A3240-05						



INTERFACE DRAC-HF

