

Technical Publication
SC-1005R0

Schematics

HF Series Generators

REVISION HISTORY

REVISION	DATE	REASON FOR CHANGE
<i>Refer to each schematic</i>		

This Document is the english original version, edited and supplied by the manufacturer.

The Revision state of this Document is indicated in the code number shown at the bottom of this page.

ADVISORY SYMBOLS

The following advisory symbols will be used throughout this manual. Their application and meaning are described below.



DANGERS ADVISE OF CONDITIONS OR SITUATIONS THAT IF NOT HEDED OR AVOIDED WILL CAUSE SERIOUS PERSONAL INJURY OR DEATH.



ADVISE OF CONDITIONS OR SITUATIONS THAT IF NOT HEDED OR AVOIDED COULD CAUSE SERIOUS PERSONAL INJURY, OR CATASTROPHIC DAMAGE OF EQUIPMENT OR DATA.




Advise of conditions or situations that if not heeded or avoided could cause personal injury or damage to equipment or data.

Note 

Alert readers to pertinent facts and conditions. Notes represent information that is important to know but which do not necessarily relate to possible injury or damage to equipment.

TABLE OF CONTENTS

SCH. No.	SCHEMATIC	REMARKS
54302010	Compact Generator	
A3274-01	Delayed Switch-Off PCB	<i>Only used with High Speed Generator</i>
A3000-xx	HT Controller PCB	<i>See Board code in Generator</i>
A3004-xx	Filament Control PCB	<i>See Board code in Generator</i>
A3009-xx	Interface Control PCB	<i>See Board code in Generator</i>
A3063-06	IPM Driver PCB	
A3212-01	Charge / Discharge Monitor PCB	
A3024-xx	ATP Console CPU PCB	
A3012-xx	AEC Control PCB	<i>Only used with AEC Option - See Board code in Generator</i>
A3129-21	RAD Display PCB	<i>RAD Console</i>




THIS SHEET ONLY APPLIES TO THREE PHASE GENERATORS MANUFACTURED FOR INPUT LINES FROM 400 VAC

Note.- In Three Phase Generators manufactured for Input Lines up to 240 VAC, Input Transformer 6T2 is as shown on sheet 2.

Note.- 12TS4, 12F14, 12F15, and 12F16 are optional

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	54302010											
H	NC 12/466	R. Asenjo	08/11/12															
G	CN 10/533	F. Garcia	14/04/11	DRAWING	F.GARCIA	30/09/02	1/8											
F	CN 10/226	F. Garcia	17/05/10	REVISED	A.DIAZ	11/11/02												
E	CN 09/385	F. Garcia	30/12/09								H	G	F	E	D	C	I	←REV
I	NC 13/015	R. Asenjo	21/01/13								SEDECAL							
							COMPACT-LV DRAC/LF RAC GENERATOR											




**THIS SHEET ONLY APPLIES TO SINGLE PHASE GENERATORS
MANUFACTURED FOR INPUT LINES UP TO 240 VAC**

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	54302010					
H	NC 12/466	R. Asenjo	08/11/12				H G F E D C I ←REV					
G	CN 10/533	F. Garcia	14/04/11	DRAWING	F.GARCIA	30/09/02	2/8					
F	CN 10/226	F. Garcia	17/05/10	REVISED	A.DIAZ	11/11/02						
E	CN 09/385	F. Garcia	30/12/09									
I	NC 13/015	R. Asenjo	21/01/13									

**COMPACT-LV DRAC/LF RAC
GENERATOR**

THIS SHEET ONLY APPLIES TO GENERATORS WITH HT TRANSFORMERS FOR TWO X-RAY TUBES




NOTE 2.- For High Speed Version:
 -T2 SCL signal is shortcircuited with -TUBE 2 SCL signal.
 -ACC and -BRAKE signals for Low Speed Version becomes -TUBE 1 RTR and -TUBE 2 RTR respectively.
 6J1 connector used only for High Speed version (see LV-DRAC documentation).

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF
H	NC 12/466	R. Asenjo	08/11/12			
G	CN 10/533	F. Garcia	14/04/11	DRAWING	F.GARCIA	30/09/02
F	CN 10/226	F. Garcia	17/05/10	REVISED	A.DIAZ	11/11/02
E	CN 09/385	F. Garcia	30/12/09			
I	NC 13/015	R. Asenjo	21/01/13			

54302010

COMPACT-LV DRAC/LF RAC GENERATOR

THIS SHEET ONLY APPLIES TO GENERATORS WITH HT TRANSFORMER FOR ONLY ONE X-RAY TUBE









NOTE 2.- For High Speed Version:
 -T2 SCL signal is shortcircuited with -TUBE 2 SCL signal.
 -ACC and -BRAKE signals for Low Speed Version becomes -TUBE 1 RTR and -TUBE 2 RTR respectively.
 6J1 connector used only for High Speed version (see LV-DRAC documentation).

NOTE 1.- Only for one FL tube:
 Wires FIL 1 RTN (P4-15) and FIL 2 RTN (P4-16), are factory inverted on the INTERFACE CONTROL Board, and 6J3-13 is disconnected and insulated.

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF
H	NC 12/466	R. Asenjo	08/11/12			
G	CN 10/533	F. Garcia	14/04/11	DRAWING	F.GARCIA	30/09/02
F	CN 10/226	F. Garcia	17/05/10	REVISED	A.DIAZ	11/11/02
E	CN 09/385	F. Garcia	30/12/09			
I	NC 13/015	R. Asenjo	21/01/13			



54302010

COMPACT-LV DRAC/LF RAC GENERATOR




REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	54302010							
H	NC 12/466	R. Asenjo	08/11/12				H	G	F	E	D	C	I	← REV
G	CN 10/533	F. Garcia	14/04/11	DRAWING	F.GARCIA	30/09/02								
F	CN 10/226	F. Garcia	17/05/10	REVISED	A.DIAZ	11/11/02								
E	CN 09/385	F. Garcia	30/12/09											
I	NC 13/015	R. Asenjo	21/01/13											



COMPACT-LV DRAC/LF RAC GENERATOR



Cable to High Speed Starter connection 6J1 (sheets 3 & 4)



Note:
Connect FAN to 11TS2-17 or 11TS2-18 as required.



THIS SHEET ONLY APPLIES TO HIGH SPEED STARTER FOR ONLY ONE X-RAY TUBE (OPTIONAL)

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	54302010					
H	NC 12/466	R. Asenjo	08/11/12				← REV					
G	CN 10/533	F. Garcia	14/04/11	DRAWING	F.GARCIA	30/09/02	6/8					
F	CN 10/226	F. Garcia	17/05/10	REVISED	A.DIAZ	11/11/02						
E	CN 09/385	F. Garcia	30/12/09									
I	NC 13/015	R. Asenjo	21/01/13									







COMPACT-LV DRAC/LF RAC GENERATOR

THIS SHEET ONLY APPLIES TO LOW SPEED STARTER

NOTE.- 10T3, 10C5-1, and 10R4-1 for 330 VAC Stator

Note.- Capacitor 10C5 could be 40 uF depending on tube used




Note.- 10J7 and 10J8 are aerial connectors



REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	54302010							
H	NC 12/466	R. Asenjo	08/11/12				H	G	F	E	D	C	I	←REV
G	CN 10/533	F. Garcia	14/04/11	DRAWING	F.GARCIA	30/09/02								
F	CN 10/226	F. Garcia	17/05/10	REVISED	A.DIAZ	11/11/02								
E	CN 09/385	F. Garcia	30/12/09											
I	NC 13/015	R. Asenjo	21/01/13											




COMPACT-LV DRAC/LF RAC GENERATOR



061109

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	DWG:
				F. GARCIA	05/05/98	1 / 1	A3274-01
				A. DIAZ	05/05/98		
C	CN 00/183	F. GARCIA	06/09/00				
B	CN 99/003	F. GARCIA	13/01/99				
							C B ← REV


DELAYED SWITCH-OFF



(* Use jumpers only for A3000-34 and 35

L	NC 09/266	F. Díaz	13/10/09	NAME	DATE	SHEET / OF	A3000-33/34/35/36/43/44/51								
K	NC 05/061	F. Díaz	01/04/05	DRAWING	F. Díaz	04/11/09	1/3	N	M	L	K	J	I	O	← REV
O	NC 12/242	F. Díaz	04/06/12	REVISED	A. Díaz	04/11/09									
N	NC 10/090	F. Díaz	24/02/10												
M	NC 09/300	F. Díaz	04/11/09												
REV	DESCRIPTION	ISSUED BY	DATE	SEDECAL			HT CONTROLLER								

HV TANK	JP3	JP4	JP5	JP6
COMPACT	2	1	2	2
VERTICAL	1	2	1	1



- (1) The value of R28 is 86K6 for A3000-35
- (2) (REV)1 - Remove D19 for three phase generators
- (3) R16, R17 = 4K02, and R20, R21 = 69K8 for A3000-36
- (7) Adjust R29 to have a period of 48 us at U19-5 for A3000-43
- (9) C41 = 100nF. Adjust R29 to have a period of 48 us at U19-5 for A3000-51

L	NC 09/266	F. Diaz	13/10/09	NAME	DATE	SHEET / OF	A3000-33/34/35/36/43/44/51	
K	NC 05/061	F. Diaz	01/04/05	DRAWING	F. Diaz	04/11/09		2/3
				REVISED	A. Diaz	04/11/09		
O	NC 12/242	F. Diaz	04/06/12					
N	NC 10/090	F. Diaz	24/02/10					
M	NC 09/300	F. Diaz	04/11/09					
REV	DESCRIPTION	ISSUED BY	DATE	SEDECAL			HT CONTROLLER	


(4) R38, R39 = 226K, and R60, R62 = 25K5 for A3000-36

(5) R36, R37 = 200K, and R40 = 2K21 for A3000-44

(8) Remove Z6 and Z7 for A3000-44

1V = 10mA (FROM 10mA UP TO 80 mA)
 1V = 100mA (FROM 100mA UP TO maximum mA)
 1V = 20mA (FROM 10mA UP TO 80 mA)
 1V = 200mA (FROM 100mA UP TO maximum mA)


A3000-44



NOTE 1: W1-A for version A3004-10/11/12, W1-B for version A3004-09
 NOTE 1: W1-A para versión A3004-10/11/12, W1-B para versión A3004-09

NOTE 2: R7= 2K67 ohms for A3004-09/10
 R7= 1K4 ohms for A3004-11
 R7= 1K ohms for A3004-12

NOTE 2: R7= 2K67 ohms para A3004-09/10
 R7= 1K4 ohms para A3004-11
 R7= 1K ohms para A3004-12




PCB 90068-04

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3004-09/10/11/12									
F	NC 09/347	F. Díaz	26/11/09	DRAWING	F. Díaz	26/10/04	1/1	F	E	Db1	Cb1	Bb1	Ab1	b2	b1	← REV
E	NC 09/162	R. Asenjo	05/06/09	REVISED	A. Díaz	26/10/04										
Db1	NC 06/251	F. Díaz	01/12/06													
Cb1	NC 05/112															



Filament Driver HC



VERSION	USED FOR	COMPONENTS ON SHEET 2	L1, L2, L3, K6
A3009-09	RAD + HV	REQUIRED	NOT REQUIRED
A3009-10	R&F	NOT REQUIRED	REQUIRED
A3009-11	ONLY RAD	NOT REQUIRED	NOT REQUIRED
A3009-12	R&F + HV	REQUIRED	REQUIRED
A3009-15	MAMO	NOT REQUIRED	NOT REQUIRED

- NOTES:
- 1.- FILAMENTO 1 IS SMALL FOR RAD AND LARGE FOR SF
 - 2.- HIGH VOLTAGE SUPPLY ON AEC OR ABS UNITS ONLY
 - 3.- L1, L2, L3 AND K6 NOT USED ON RAD SYSTEM
 - 4.- K5 AND D11 NOT PRESENT IN A3009-15
- NOTAS:
- 1.- FILAMENTO 1 ES FINO PARA RAD Y GRUESO PARA SF
 - 2.- LA FUENTE DE ALTA TENSION SOLO PARA UNIDADES CON AEC O ABS
 - 3.- L1, L2, L3 Y K6 NO SE USAN EN SISTEMAS RAD
 - 4.- K5 Y D11 NO SE MONTAN EN A3009-15

REV	DESCRIPTION	ISSUED BY	DATE	DRAWING	NAME	DATE	SHEET / OF	A3009-09/10/11/12/15						
				REVISED	A. Díaz	14/02/06	1/2	M	L	K	Jb1	I b2	I b1	← REV
M	NC 13/098	M.González	13/03/13											
L	NC 09/360	F. Díaz	30/11/09											
K	NC 08/102	F. Díaz	23/04/08											
								INTERFACE CONTROL						




PHOTO TUBE SPLY OUTPUT (P2-1)

+500 VDC OR +300 VDC : JUMPERS W3 TO W8 IN POS "A", AND ADJUSTING WITH R20
 VARIABLE NEGATIVE OUTPUT : JUMPERS W3 TO W10 IN POS "B"


SALIDA PHOTO TUBE SPLY (P2-1)

+500 VDC O +300 VDC : PUENTES W3 A W8 EN POS "A",Y AJUSTAR CON R20
 SALIDA NEGATIVA VARIABLE : PUENTES W3 A W10 EN POS "B"

NOTE.- BE SURE THAT THE VALUE OF THE PHOTO TUBE SPLY OUTPUT IN P2-1 IS THE REQUIRED BEFORE CONECTING THE SYSTEM INTERFACE

NOTA.- ASEGURESE DE QUE EL VALOR DE LA SALIDA PHOTO TUBE SPLY EN P2-1 ES EL REQUERIDO ANTES DE HACER LA INTERCONEXION DEL SISTEMA


REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3009-09/10/11/12/15								
				DRAWING	F. Díaz	14/02/06	2/2		M	L	K	Jb1	I b2	I b1	← REV
				REVISED	A. Díaz	14/02/06									
M	NC 13/098	M.González	13/03/13												
L	NC 09/360	F. Díaz	30/11/09												
K	NC 08/102	F. Díaz	23/04/08												
SEDECAL												INTERFACE CONTROL			




C.I. 90132-06

				NAME	DATE	SHEET / OF	A3063-06		
				DRAWING	F. Diaz	20/06/11	1/1		
				REVISED	A. Diaz	20/06/11			
							IPM DRIVER		
B	NC 13/001	F. Diaz	04/01/13						
REV	DESCRIPTION	ISSUED BY	DATE						

TO MAIN STORAGE CAPACITORS




E	CN 98/114	F.GARCIA	10/10/98	CHARGE/DISCHARGE MONITOR SEDECAL S.A.	ENG	F. GARCIA
D	CN 98/057	F.GARCIA	29/04/98		CHK	A. DIAZ
C	CN 96/033	F.GARCIA	06/03/96		REV	CN99/80 (22/07/99)
B	CN 95/080	F.GARCIA	10/05/95		DWG	3212-01 REV. F
A	CN 95/74	F.GARCIA	25/04/95			REV.
REV	DESCRIPTION	ISSUED BY	DATE			REV.
				DATE		11/01/95



S	NC 12/278	F. Diaz	25/06/12	NAME	DATE	SHEET / OF	A-3024-XX						
R	NC 11/427	J.A.Garcia	25/11/11	DRAWING	F. Garcia	26/03/09	1/4						
Q	NC 11/036	F. Diaz	09/02/11	REVISED	A. Diaz	26/03/09							
P	add version -97	M.A. González	25/05/10					S	R	Q	P	O	← REV
P	NC 08/269	M.A. González	11/09/08										
O	NC 07/103	F. Diaz	07/05/07										ATP CONSOLE
REV	DESCRIPTION	ISSUED BY	DATE										

NOTE.- ALL DIODES 1N4148 UNLESS SPECIFIED






JUMPER	NORMAL	MOBILE
JP14	set	removed
JP15	removed	set
JP16	removed	set
JP17	removed	set
JP18	removed	set

note: only for -97 version


S	NC 12/278	F. Diaz	25/06/12	NAME	DATE	SHEET / OF	A-3024-XX						
R	NC 11/427	J.A.Garcia	25/11/11	DRAWING	F. Garcia	26/03/09	3/4	S	R	Q	P	O	← REV
Q	NC 11/036	F. Diaz	09/02/11	REVISED	A. Diaz	26/03/09							
P	add version -97	M.A. González	25/05/10										
P	NC 08/269	M.A. González	11/09/08										
O	NC 07/103	F. Diaz	07/05/07										
REV	DESCRIPTION	ISSUED BY	DATE	SEDECAL		ATP CONSOLE							

SERIAL COMMUNICATION REQUIRES THE COMPONENTS SHOWN IN THIS BOX

version -26: RS422 communication, JP7-JP8-JP9-JP10-JP11-JP21-JP22 "B" position
 version -27: RS485 communication, JP7-JP8-JP21-JP22 "A" position; JP9, JP10, JP11 "B" position.





S	NC 12/278	F. Diaz	25/06/12	NAME	DATE	SHEET / OF	A-3024-XX						
R	NC 11/427	J.A.Garcia	25/11/11	DRAWING	F. Garcia	26/03/09	4/4	S	R	Q	P	O	← REV
Q	NC 11/036	F. Diaz	09/02/11	REVISED	A. Diaz	26/03/09							
P	add version -97	M.A. González	25/05/10										
P	NC 08/269	M.A. González	11/09/08										
O	NC 07/103	F. Diaz	07/05/07										
REV	DESCRIPTION	ISSUED BY	DATE	SEDECAL		ATP CONSOLE							



BOARD	VERSION	JP1	JP2	JP4
3012-01	PHOTOMULTIPLIER	B	B	A
3012-02	FOR TV CAMERA	A	A	A
3012-05	4 ION CHAMBERS	C	A	B


NOTE.- JUMPER JP3 IN POSITION A FOR ION CHAMBERS WITH HIGH SENSIBILITY (GAIN 2V / mR)

AEC CONTROL SEDECAL S.A.		ENG	F. GARCIA
		CHK	A. DIAZ
DWG	REV	CN 08/137 (14/05/08)	
		A3012-01/02	REV. L
		A3012-05	REV. E
DATE	15/05/08		




R5, R6, R7 and R8 are 82 ohm, 2w
 R24, R25, R26 and R27 are 47 ohm, 2w


NAME		DATE		SHEET / OF		A3129-21	
DRAWING		F.GARCIA		05/05/99		1/3	
REVISED		A.DIAZ		12/12/99		A ← REV	
A		NC: 12/210		J.A.Garcia		25/05/12	
REV	DESCRIPTION	ISSUED BY	DATE			RAD DISPLAY	



NOTE.- SEVEN SEGMENT DISPLAYS TYPE HDSP3600



		NAME	DATE	SHEET / OF	A3129-21	
		DRAWING	F.GARCIA	05/05/99	2/3	
		REVISED	A.DIAZ	12/12/99	A ← REV	
				RAD DISPLAY		
A	NC: 12/210	J.A.Garcia	25/05/12			
REV	DESCRIPTION	ISSUED BY	DATE			



NAME		DATE		SHEET / OF		A3129-21	
DRAWING		F.GARCIA		05/05/99		3/3	
REVISED		A.DIAZ		12/12/99		A ← REV	
REV		DESCRIPTION		ISSUED BY		DATE	
A	NC: 12/210	J.A.Garcia	25/05/12				



RAD DISPLAY