

Technical Publication
SC-1006R0

Schematics

HF Series Generators
Battery Generators

REVISION HISTORY

REVISION	DATE	REASON FOR CHANGE
<i>Refer to each schematic</i>		

This Document is the english original version, edited and supplied by the manufacturer.

The Revision state of this Document is indicated in the code number shown at the bottom of this page.

ADVISORY SYMBOLS

The following advisory symbols will be used throughout this manual. Their application and meaning are described below.



DANGERS ADVISE OF CONDITIONS OR SITUATIONS THAT IF NOT HEEDED OR AVOIDED WILL CAUSE SERIOUS PERSONAL INJURY OR DEATH.



ADVISE OF CONDITIONS OR SITUATIONS THAT IF NOT HEEDED OR AVOIDED COULD CAUSE SERIOUS PERSONAL INJURY, OR CATASTROPHIC DAMAGE OF EQUIPMENT OR DATA.



Advise of conditions or situations that if not heeded or avoided could cause personal injury or damage to equipment or data.


Note 

Alert readers to pertinent facts and conditions. Notes represent information that is important to know but which do not necessarily relate to possible injury or damage to equipment.

TABLE OF CONTENTS

SCH. No.	SCHEMATIC	REMARKS
54302013	ESM Compact Generator w/o Stand-Alone	<i>Generator without Stand-Alone Option</i>
A3096-02	LF-RAC PCB	<i>Only used with Low Speed Generator</i>
A3000-xx	HT Controller PCB	<i>See Board code in Generator</i>
A3004-xx	Filament Control PCB	<i>See Board code in Generator</i>
A3009-xx	Interface Control PCB	<i>See Board code in Generator</i>
A3063-06	IPM Driver PCB	
A3212-01	Charge / Discharge Monitor PCB	
A3024-xx	ATP Console CPU PCB	
A3012-xx	AEC Control PCB	<i>Only used with AEC Option - See Board code in Generator</i>
ESM (ENERGY STORAGE MODULE) WITHOUT STAND-ALONE OPTION		
A3285-02	Battery Charger PCB	
A3580-01	Battery Monitor PCB	
A3139-01	Line Monitor PCB	
CONTROL CONSOLE		
A3553-02	Serial Console Control PCB	
A3585-21	RAD Multiuse H2 Console Display PCB	


SERIAL OPERATOR CONSOLE



NOTES:
 RXD ON GENERATOR IS CONNECTED TO TXD ON SERIAL OPERATOR CONSOLE
 TXD ON GENERATOR IS CONNECTED TO RXD ON SERIAL OPERATOR CONSOLE

GENERATOR CABINET


NOTE: FOR GENERAL INTERCONNECTIONS REFER TO SCHEMATIC 54301045 (COMPACT / COMPACT-ESM GENERATORS)



NOTE: REFER TO SCHEMATIC A6188-02 FOR RS-232/422/485 SERIAL COMMUNICATION

SERIAL INTERCONNECTION CABLE FOR SERIAL CONSOLE (A7066-xx or A3149-01)

TPC or PC INTERFACE BOX




NOTES:
 RXD AND TXD ARE INTERNALLY REVERSED IN TPC OR PC INTERFACE BOX
 REFER TO SCHEMATICS I/F-036 FOR OTHER PC INTERFACE BOX CONNECTIONS

SERIAL INTERCONNECTION CABLE FOR TPC (A7067-xx or A3352-01)
 SERIAL INTERCONNECTION CABLE FOR INTERFACE BOX (A3352-01)

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	DWG:
				F. GARCIA	07/07/99	1 / 1	54301052
				A. DIAZ	24/01/00		C B A ← REV
C	NC 03 / 050	F. GARCIA	08/03/03				
B	New schematic	F. GARCIA	08/01/02				
A	Connections	F. GARCIA	02/02/01				
<p style="text-align: center;">SERIAL COMMUNICATION TO GENERATOR SYSTEM INTERCONNECTION</p>							

ATP CONSOLE BD. (Multilayer)



REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	DWG:
				F. GARCIA	04/04/04	1 / 1	A6188-03
				A. DIAZ	04/04/04		
A	CN 04/148	F. GARCIA	09/09/04				RS-232/422/485 SERIAL COMMUNICATION



FLURO CPU BOARD

CONNECTOR J4		
SIGNAL	I/O	PIN
C-FL DAT	I	1
GND	I/O	2
-CAM FL EXP	O	3
C-FL CLK	I	4
FL-C DAT	I	5
GND	I/O	6
FL-C CLK	O	7
CAM SYNC	O	8
+12V ISO	I	9
N.U.		10

CONNECTOR J1 / J10		
SIGNAL	I/O	PIN
-4 IN SEL	O	1
-9 IN SEL	O	2
+12 VDC	O	3
V SYNC	I	4
-FT SW CMD	I	5
FL START	O	6
-CAM FL EXP	O	7
N.U.		8
-ABS	O	9
-6 IN SEL	O	10
BEEP	O	11
GND	O	12
-MEM EN	O	13
-MEM GATE	O	14
N.U.		15

ATP CONSOLE BOARD

CONNECTOR J4		
SIGNAL	I/O	PIN
C-FL DAT	O	1
-CAM FL EXP	I	2
FL-C DAT	I	3
FL-C CLK	I	4
+12V ISO	I	5
GND	I/O	6
C-FL CLK	I	7
GND	I/O	8
CAM SYNC	O	9

CONNECTOR J5		
SIGNAL	I/O	PIN
IC1 INPUT	I	1
IC3 INPUT	I	2
-FLD1 DR	O	3
-FLD3 DR	O	4
GND	O	5
IC2 INPUT	I	6
IC4 INPUT	I	7
-FLD2 DR	O	8
-STRT DR	O	9

CONNECTOR J13		
SIGNAL	I/O	PIN
TIME1	I	1
TIME2	I	2
TIME3	I	3
TIME4	I	4
-PS (DSI) SEL	I	5
-CINE (DSA) SEL	I	6
-HCF SEL	I	7
SPARE IN	I	8
-TOMO PREP	O	9
-TOMO EXP	O	10
TOMO ON	O	11
EXT REF	I	12
GND	I/O	13
EXP STOP	I	14
-FS	O	15

CONNECTOR J2		
SIGNAL	I/O	PIN
-GEN OK	O	1
-SFD SEL	O	2
-THERMOSTAT 1	I	3
-SF PREP	I	4
GND (THERM. COM)	I/O	5
-COLLIMATOR	I	6
TABLE ERR (COMP)	I	7
-ROOM LIGHT	O	8
-READY	O	9
EXP OK	O	10
ABC OUT - LEFT	O	11
-DIRECT SEL	O	12
PT INPUT	I	13
EXT SYNC (FL DSI)	I	14
SPARE IN2	I	15
-SFC (-PT SEL)	I	16
-FL EXP	I	17
GND (DOOR RTN)	I/O	18
-DOOR	I	19
-THERMOSTAT 2	I	20
-AUTO OFF	O	21
SPARE IN1	I	22
ALOE	O	23
-ALOE	O	24
-ACT EXP	O	25

CONNECTOR J1		
SIGNAL	I/O	PIN
GND UNR	I	1
+12V UNR	I	2
+12VDC	O	3
-PWR OFF	O	4
-PWR ON	O	5
GND UNR	I	6
GND	I	7
+12V UNR	I	8
-12VDC	I	9

CONNECTOR J3		
SIGNAL	I/O	PIN
-EXP	O	1
C-HT DAT	O	2
HT-C DAT	I	3
-PREP	O	4
C-HT CLK	O	5
HT-C CLK	I	6
HV PT CRL	O	7
-KV DWN	O	8
LINE SYNC	I	9
-LINE CONT	O	10
-BUCKY 2 DR CMD	O	11
-BUCKY EXP	I	12
-HT INL	I	13
-BUCKY 1 DR CMD	O	14
-KV UP	O	15

HAND-SWITCH (or VET PEDAL SWITCH)

CONNECTOR J15		
SIGNAL	I/O	PIN
COM	O	1
PREP	O	2
EXP	I	3
N.U.		4

HAND-SWITCH (or VET PEDAL SWITCH)



FLURO CABLE

AEC CABLE

TOMO CABLE

INTERFACE CABLE

POWER CABLE

COMMUNICATION CABLE

GND CABLE

ADAPTATIONS BOARDS



NOTE :
Signal for Thermostats go to 4TS3 and then to TS2 or go directly to TS2 depending on Generator model

TERMINAL BLOCK 3TS1		
PIN	I/O	SIGNAL
1	O	SUPPLY (BUCKY 1)
2	O	BUCKY SPLY 1
3	O	0 VAC (BUCKY 1)
4	O	BUCKY 1 DR
5	I	-BUCKY 1 MOT
6	I	BUCKY 1 MOT RTN
7	O	GND
8	O	SUPPLY (BUCKY 2)
9	O	BUCKY SPLY 2
10	O	0 VAC (BUCKY 2)
11	O	BUCKY 2 DR
12	I	-BUCKY 2 MOT
13	I	BUCKY 2 MOT RTN
14	O	GND
15	O	+12V
16	O	-12V
17	O	+24V UNR
18	O	GND
22	I	-DOOR
23	O	DOOR RTN
24	O	-ROOM LIGHT
26	O	220 VAC SW
27	O	115 VAC SW
36	O	-SF PREP
37	O	-FL EXP
39	O	PT SPLY
42	O	IC GND
47	O	ROOM LIGHT SUP
48	O	ROOM LIGHT SW
51	I	PT INPUT
52	O	-ALOE
53	O	-READY
54	O	220 VAC SPLY

TERMINAL BLOCK 4TS3		
PIN	I/O	SIGNAL
1	I	-THERMOSTAT 1
2	I	-THERMOSTAT 2
3	I	THERM. COMM
5	I	TABLE ERROR
7	I	-SFC (-PT SEL)
8	O	ALOE
10	O	-GEN OK
11	O	-SFD SEL
12	O	-DIRECT SEL
15	O	EXP OK
16	O	-ACT EXP
17	O	-AUTO OFF
18	I	ABC OUT - LEFT
19	I	EXT SYNC (FL DSI)
20	I	-COLLIMATOR
A	I	Foot Switch Cmd.
B		GND
C	O	Prep/Rdy Acq Rad
D	I	Boost Fluoro
E	I	Digital Exposure

CONNECTOR 6J3		
PIN	I/O	SIGNAL
1	O	-EXP
2	O	C-HT DAT
3	I	HT-C DAT
4	O	-PREP
5	O	C-HT CLK
6	I	HT-C CLK
7	O	HV PT CRL
8	O	-KV DWN
9	I	LINE SYNC
10	O	-LINE CONT
11	O	-BUCKY 2 DR CMD
12	I	-BUCKY EXP
13	I	-HT INL
14	O	-BUCKY 1 DR CMD
15	O	-KV UP
16	O	GND

COMPACT GENERATOR CABINET



NOTE - For Generator with DRAC :
The ROTOR TUBE connections are made to TS2 on the DRAC

TERMINAL BLOCK 4TS2 / 10TS2 / 11TS2		
SIGNAL	I/O	PIN
MAIN T1	O	1
AUX T1	O	2
COM T1	O	3
-THERMOSTAT 1	I	4
THERMOST. COMM.	I	5
FAN 1	O	6
0 VAC	O	7
GND	O	8
MAIN T2	O	9
AUX T2	O	10
COM T2	O	11
-THERMOSTAT 2	I	12
THERMOST. COMM.	I	13
FAN 2	O	14
0 VAC	O	15
GND	O	16

NOTE :
TS2-6 is GND when TS2 only has 6 terminals.

CONNECTOR 6J2		
PIN	I/O	SIGNAL
1	O	GND
2	O	+12V UNR
3	O	+12VDC
4	I	-PWR OFF
5	I	-PWR ON
6	O	GND
7	O	GND
8	O	+12V UNR
9	O	-12VDC

NOTE - For Serial Cabinet :
- the Interface cable connections are factory made to 3TS1, TS2 and 4TS3 terminal blocks.
- 6J2 is not supplied. Power cable is directly connected in factory to J1 of the ATP Console.
- 6J3 is not supplied. Power cable is directly connected in factory to J3 of the ATP Console.

LOCKS BOARD TERMINAL BLOCK TB7		
SIGNAL	I/O	PIN
24 VAC LAMP	O	3
0 VAC LAMP	O	4
+24 VDC LOCKS	O	5
0 VDC LOCKS	O	6
+24 VDC LOCKS	O	7
0 VDC LOCKS	O	8

COLLIMATOR LAMP
LOCKS (Table & Tube Stand)


REV	DESCRIPTION	ISSUED BY	DATE
F	Adaptation Boards	F. GARCIA	01/09/06
E	CN 04/032	F. GARCIA	26/04/04
D	New interface	F. GARCIA	08/01/02
C	New interface	F. GARCIA	02/02/01
B	CN 00/211	F. GARCIA	04/01/01
A	New interface	F. GARCIA	05/05/00

NAME	DATE	SHEET / OF
DRAWING	F. GARCIA	04/04/99
REVISED	A. DIAZ	28/01/00

DWG:	54301045			
	F	E	D	C
				← REV




COMPACT / COMPACT-ESM GENERATORS SYSTEM INTERCONNECTION



REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF
I	NC 13/015	R. Asenjo	21/01/13			
H	NC 12/466	R. Asenjo	08/11/12	DRAWING	F.GARCIA	14/03/05
G	CN 10/533	F. DIAZ	14/04/11	REVISED	A.DIAZ	14/03/05
F	CN 10/226	F. GARCIA	17/05/10			
E	CN 06/159	F. GARCIA	08/08/06			
D	NE 06/166	F. GARCIA	07/07/06			

SEDECAL						
54302013						
PSU (without Stand-Alone)						
ESM COMPACT-LV DRAC/LF RAC GENERATOR						
(ONE/TWO X-RAY TUBES)						

Note:
See DELAYED SWITCH-OFF
board on sheet 5 for the
connections (*) added to J4



Cable to High Speed Starter connection 6J1 (sheets 2 & 3)


**CONTROL DRAC
PCB (A3243-03)**

11A1

See Note

J4

- 1 +5 VDC
- 2 +5 VDC
- 3 +5 VDC
- 4 +5 VDC
- 5 GND
- 6 GND
- 7 GND
- 8 GND
- 9 +15 VDC
- 10 +15 VDC
- 11 +15 VDC
- 12 +15 VDC
- 13 -15 VDC
- 14 -15 VDC
- 15 V UNR
- 16 V UNR
- 17 PRECH IF
- 18 PRECH IF
- 19 +5 VDC
- 20 +5 VDC
- 21 T3 IF
- 22 T3 IF
- 23 T2 IF
- 24 T2 IF
- 25 T1 IF
- 26 T1 IF




Ribbon Cable



**INTERFASE DRAC-HF
PCB (A3240-05)**

11A2


- 1 +5 VDC
- 2 +5 VDC
- 3 +5 VDC
- 4 +5 VDC
- 5 GND
- 6 GND
- 7 GND
- 8 GND
- 9 +15 VDC
- 10 +15 VDC
- 11 +15 VDC
- 12 +15 VDC
- 13 -15 VDC
- 14 -15 VDC
- 15 V UNR
- 16 V UNR
- 17 PRECH IF
- 18 PRECH IF
- 19 +5 VDC
- 20 +5 VDC
- 21 T3 IF
- 22 T3 IF
- 23 T2 IF
- 24 T2 IF
- 25 T1 IF
- 26 T1 IF



Note:
Connect FAN to 11TS2-17
or 11TS2-18 as required.



Note.- 11J7 is an aerial connector




THIS SHEET ONLY APPLIES TO HIGH SPEED STARTER FOR ONLY ONE X-RAY TUBE (OPTIONAL)

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	54302013						
I	NC 13/015	R. Asenjo	21/01/13				5/9						
H	NC 12/466	R. Asenjo	08/11/12	DRAWING	F.GARCIA	14/03/05	← REV						
G	CN 10/533	F. DIAZ	14/04/11	REVISED	A.DIAZ	14/03/05	H	G	F	E	D	C	I
F	CN 10/226	F. GARCIA	17/05/10										
E	CN 06/159	F. GARCIA	08/08/06										
D	NE 06/166	F. GARCIA	07/07/06										






PSU (without Stand-Alone)
ESM COMPACT-LV DRAC/LF RAC GENERATOR
(ONE/TWO X-RAY TUBES)



Cable to High Speed Starter connection 6J1 (sheets 3 & 4)

Note.- 11J7 is an aerial connector



THIS SHEET ONLY APPLIES TO HIGH SPEED STARTER FOR ONLY TWO X-RAY TUBES (STANDARD)

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF
I	NC 13/015	R. Asenjo	21/01/13			
H	NC 12/466	R. Asenjo	08/11/12	DRAWING	F.GARCIA	14/03/05
G	CN 10/533	F. DIAZ	14/04/11	REVISED	A.DIAZ	14/03/05
F	CN 10/226	F. GARCIA	17/05/10			6/9
E	CN 06/159	F. GARCIA	08/08/06			
D	NE 06/166	F. GARCIA	07/07/06			

54302013						
H	G	F	E	D	C	I ← REV



SEDECAL	
PSU (without Stand-Alone) ESM COMPACT-LV DRAC/LF RAC GENERATOR (ONE/TWO X-RAY TUBES)	

THIS SHEET ONLY APPLIES TO LOW SPEED STARTER






NOTE.- 10T3, 10C5-1, and 10R4-1 for 330 VAC Stator

Note.- Capacitor 10C5 could be 40 uF depending on tube used

Note.- 10J7 and 10J8 are aerial connectors




NOTE.- Use Tube 2 connections for FL tube





REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF
I	NC 13/015	R. Asenjo	21/01/13			
H	NC 12/466	R. Asenjo	08/11/12	DRAWING	F.GARCIA	14/03/05
G	CN 10/533	F. DIAZ	14/04/11	REVISED	A.DIAZ	14/03/05
F	CN 10/226	F. GARCIA	17/05/10			
E	CN 06/159	F. GARCIA	08/08/06			
D	NE 06/166	F. GARCIA	07/07/06			

54302013 PSU (without Stand-Alone) ESM COMPACT-LV DRAC/LF RAC GENERATOR (ONE/TWO X-RAY TUBES)						







NOTE:
 9 Amp. BATTERIES




0 VAC J2-5
 20 VAC J2-6
 Input Line J2-3

to J2 connector of the LINE MONITOR board (sheet 4)

* Connect according to the Input Voltage Line



REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	54302013					
I	NC 13/015	R. Asenjo	21/01/13									
H	NC 12/466	R. Asenjo	08/11/12	DRAWING	F.GARCIA	14/03/05						
G	CN 10/533	F. DIAZ	14/04/11	REVISED	A.DIAZ	14/03/05						
F	CN 10/226	F. GARCIA	17/05/10									
E	CN 06/159	F. GARCIA	08/08/06									
D	NE 06/166	F. GARCIA	07/07/06									
							SEDECAL			PSU (without Stand-Alone) ESM COMPACT-LV DRAC/LF RAC GENERATOR (ONE/TWO X-RAY TUBES)		



NOTE:
17 Amp. BATTERIES

**BATTERY CHARGER BOARD
(A3285-02)**

**BAT. CHARGER TRANSFORMER
2T1**

**BATTERY MONITOR BOARD
(A3580-01)**


REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF
I	NC 13/015	R. Asenjo	21/01/13			
H	NC 12/466	R. Asenjo	08/11/12	DRAWING	F.GARCIA	14/03/05
G	CN 10/533	F. DIAZ	14/04/11	REVISED	A.DIAZ	14/03/05
F	CN 10/226	F. GARCIA	17/05/10			
E	CN 06/159	F. GARCIA	08/08/06			
D	NE 06/166	F. GARCIA	07/07/06			



54302013


H G F E D C I ←REV


PSU (without Stand-Alone)
ESM COMPACT-LV DRAC/LF RAC GENERATOR
(ONE/TWO X-RAY TUBES)



Note:
 The LF-RAC shown in the Generator schematic (543020XX) is a basic documentation, only for interconnection purposes.
 Follow this schematic for detailed component information of the LF-RAC board.

061109

				NAME	DATE	SHEET / OF	DWG: A3096-02			
				DRAWING	F. GARCIA	21/09/97				
				REVISED	A. DIAZ	21/09/97				
							1 / 1			
							B ← REV			
							LF - RAC Board (Line Frequency Rotatory Anode Controller)			
B	CN 03/004	F. Garcia	13/01/03							
REV	DESCRIPTION	ISSUED BY	DATE							



(* Use jumpers only for A3000-34 and 35


(*

(*

(*

L	NC 09/266	F. Diaz	13/10/09	NAME	DATE	SHEET / OF	A3000-33/34/35/36/43/44/51								
P	NC13/315	R. Herosilla	30/09/13	DRAWING	F. Diaz	04/11/09	1/3	N	M	L	K	J	P	O	← REV
O	NC 12/242	F. Diaz	04/06/12	REVISED	A. Diaz	04/11/09									
N	NC 10/090	F. Diaz	24/02/10												
M	NC 09/300	F. Diaz	04/11/09												
REV	DESCRIPTION	ISSUED BY	DATE				HT CONTROLLER								

HV TANK	JP3	JP4	JP5	JP6
COMPACT	2	1	2	2
VERTICAL	1	2	1	1



- (1) The value of R28 is 86K6 for A3000-35
- (2) (REV)1 - Remove D19 for three phase generators
- (3) R16, R17 = 4K02, and R20, R21 = 69K8 for A3000-36
- (7) Adjust R29 to have a period of 48 us at U19-5 for A3000-43
- (9) C41 = 100nF. Adjust R29 to have a period of 48 us at U19-5 for A3000-51


L	NC 09/266	F. Diaz	13/10/09	NAME	DATE	SHEET / OF	A3000-33/34/35/36/43/44/51									
P	NC13/315	R. Hermosilla	30/09/13	DRAWING	F. Diaz	04/11/09	2/3		N	M	L	K	J	P	O	← REV
O	NC 12/242	F. Diaz	04/06/12	REVISED	A. Diaz	04/11/09										
N	NC 10/090	F. Diaz	24/02/10													
M	NC 09/300	F. Diaz	04/11/09													
REV	DESCRIPTION	ISSUED BY	DATE					HT CONTROLLER								

(4) R38, R39 = 226K, and R60, R62 = 25K5 for A3000-36

(5) R36, R37 = 200K, and R40 = 2K21 for A3000-44

(8) Remove Z6 and Z7 for A3000-44

1V = 10mA (FROM 10mA UP TO 80 mA)
 1V = 100mA (FROM 100mA UP TO maximum mA)
 1V = 20mA (FROM 10mA UP TO 80 mA)
 1V = 200mA (FROM 100mA UP TO maximum mA)




L	NC 09/266	F. Diaz	13/10/09	NAME	DATE	SHEET / OF	A3000-33/34/35/36/43/44/51								
				DRAWING	F. Diaz	04/11/09	3/3	N	M	L	K	J	P	O	← REV
P	NC13/315	R. Hermosilla	30/09/13	REVISED	A. Diaz	04/11/09									
O	NC 12/242	F. Diaz	04/06/12												
N	NC 10/090	F. Diaz	24/02/10												
M	NC 09/300	F. Diaz	04/11/09												
REV	DESCRIPTION	ISSUED BY	DATE					HT CONTROLLER							

NOTE 1: W1-A for version A3004-10/11/12, W1-B for version A3004-09
 NOTE 1: W1-A para versión A3004-10/11/12, W1-B para versión A3004-09

NOTE 2: R7= 2K67 ohms for A3004-09/10
 R7= 1K4 ohms for A3004-11
 R7= 1K ohms for A3004-12

NOTE 2: R7= 2K67 ohms para A3004-09/10
 R7= 1K4 ohms para A3004-11
 R7= 1K ohms para A3004-12




PCB 90068-04

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3004-09/10/11/12									
				DRAWING	F. Díaz	26/10/04	1/1	F	E	Db1	Cb1	Bb1	Ab1	b2	b1	← REV
				REVISED	A. Díaz	26/10/04										
					NC 09/347	F. Díaz	26/11/09									
					NC 09/162	R. Asenjo	05/06/09									
					NC 06/251	F. Díaz	01/12/06									
					NC 05/112											



Filament Driver HC



VERSION	USED FOR	COMPONENTS ON SHEET 2	L1, L2, L3, K6
A3009-09	RAD + HV	REQUIRED	NOT REQUIRED
A3009-10	R&F	NOT REQUIRED	REQUIRED
A3009-11	ONLY RAD	NOT REQUIRED	NOT REQUIRED
A3009-12	R&F + HV	REQUIRED	REQUIRED
A3009-15	MAMO	NOT REQUIRED	NOT REQUIRED

NOTES:
 1.- FILAMENT 1 IS SMALL FOR RAD AND LARGE FOR SF
 2.- HIGH VOLTAGE SUPPLY ON AEC OR ABS UNITS ONLY
 3.- L1, L2, L3 AND K6 NOT USED ON RAD SYSTEM
 4.- K5 AND D11 NOT PRESENT IN A3009-15

NOTAS:
 1.- FILAMENTO 1 ES FINO PARA RAD Y GRUESO PARA SF
 2.- LA FUENTE DE ALTA TENSION SOLO PARA UNIDADES CON AEC O ABS
 3.- L1, L2, L3 Y K6 NO SE USAN EN SISTEMAS RAD
 4.- K5 Y D11 NO SE MONTAN EN A3009-15

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3009-09/10/11/12/15				
N	NC 14/118	M.González	16/04/14	DRAWING	F. Díaz	14/02/06	1/2				
M	NC 13/098	M.González	13/03/13	REVISED	A. Díaz	14/02/06					
L	NC 09/360	F. Díaz	30/11/09								
K	NC 08/102	F. Díaz	23/04/08								



INTERFACE CONTROL




PHOTO TUBE SPLY OUTPUT (P2-1)

+500 VDC OR +300 VDC : JUMPERS W3 TO W8 IN POS "A", AND ADJUSTING WITH R20
 VARIABLE NEGATIVE OUTPUT : JUMPERS W3 TO W10 IN POS "B"


SALIDA PHOTO TUBE SPLY (P2-1)

+500 VDC O +300 VDC : PUENTES W3 A W8 EN POS "A",Y AJUSTAR CON R20
 SALIDA NEGATIVA VARIABLE : PUENTES W3 A W10 EN POS "B"


NOTE.- BE SURE THAT THE VALUE OF THE PHOTO TUBE SPLY OUTPUT IN P2-1 IS THE REQUIRED BEFORE CONECTING THE SYSTEM INTERFACE
 NOTA.- ASEGURESE DE QUE EL VALOR DE LA SALIDA PHOTO TUBE SPLY EN P2-1 ES EL REQUERIDO ANTES DE HACER LA INTERCONEXION DEL SISTEMA

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3009-09/10/11/12/15								
				DRAWING	F. Díaz	14/02/06	2/2	N	M	L	K	Jb1	l b2	l b1	←REV
				REVISED	A. Díaz	14/02/06									
												INTERFACE CONTROL			


TO MAIN STORAGE CAPACITORS



E	CN 98/114	F.GARCIA	10/10/98	CHARGE/DISCHARGE MONITOR SEDECAL S.A.	ENG	F. GARCIA
D	CN 98/057	F.GARCIA	29/04/98		CHK	A. DIAZ
C	CN 96/033	F.GARCIA	06/03/96		REV	CN99/80 (22/07/99)
B	CN 95/080	F.GARCIA	10/05/95		DWG	3212-01 REV. F
A	CN 95/74	F.GARCIA	25/04/95			REV.
REV	DESCRIPTION	ISSUED BY	DATE			REV.
				DATE		11/01/95



S	NC 12/278	F. Diaz	25/06/12	NAME	DATE	SHEET / OF	A-3024-XX								
R	NC 11/427	J.A.Garcia	25/11/11	DRAWING	F. Garcia	26/03/09	1/4	U	T	S	R	Q	P	O	←REV
Q	NC 11/036	F. Diaz	09/02/11	REVISED	A. Diaz	26/03/09									
U	NC 13/257	R. Hermosilla	19/02/14												
T	NC 13/417	R. Hermosilla	21/01/14												
REV	DESCRIPTION	ISSUED BY	DATE				ATP CONSOLE								



JP1, JP2, AND JP3 ON POSITION "A" FOR NOVEL CONSOLE

-AUTO OFF to J2-21 (page 3)

J14
1 COMM
2 PREP
3 EXP
4 POWER OFF
5 POWER ON
6 V (ON)


J9
11-12 V (ON)
1-2 COL 0
3-4 COL 1
5-6 COL 2
7-8 COL 3
9-10 COL 4
13-14 COL 5
23-24 ROW 0
27-28 ROW 1
29-30 ROW 2
25-26 ROW 3
15-16 ROW 4
17-18 ROW 5
19-20 ROW 6
21-22 ROW 7
39-40 POWER OFF
37-38 POWER ON
33-34

AEC CONTROL BOARD
(A3012-01/02/05/06)

JP5 POSITION:
A - VD SELECTION
B - IC4 SELECTION
C - PT INPUT SELECTION

S	NC 12/278	F. Díaz	25/06/12	NAME	DATE	SHEET / OF	A-3024-XX						
R	NC 11/427	J.A. Garcia	25/11/11	DRAWING	F. García	26/03/09	2/4						
Q	NC 11/036	F. Díaz	09/02/11	REVISED	A. Díaz	26/03/09	U T S R Q P O ← REV						
U	NC 13/257	R. Hermosilla	19/02/14	SEDECAL				ATP CONSOLE					
T	NC 13/417	R. Hermosilla	21/01/14										
REV	DESCRIPTION	ISSUED BY	DATE										

NOTE - ALL DIODES 1N4148 UNLESS SPECIFIED




JUMPER	NORMAL	MOBILE
JP14	set	removed
JP15	removed	set
JP16	removed	set
JP17	removed	set
JP18	removed	set

note: only for -97 version

S	NC 12/278	F. Díaz	25/06/12	NAME	DATE	SHEET / OF	A-3024-XX								
R	NC 11/427	J.A.García	25/11/11	DRAWING	F. García	26/03/09	3/4	U	T	S	R	Q	P	O	← REV
Q	NC 11/036	F. Díaz	09/02/11	REVISED	A. Díaz	26/03/09									
U	NC 13/257	R. Hermosilla	19/02/14												
T	NC 13/417	R. Hermosilla	21/01/14												
REV	DESCRIPTION	ISSUED BY	DATE				SEDECAL		ATP CONSOLE						

SERIAL COMMUNICATION REQUIRES THE COMPONENTS SHOWN IN THIS BOX



version -26: RS422 communication, JP7-JP8-JP9-JP10-JP11-JP21-JP22 "B" position
 version -27: RS485 communication, JP7-JP8-JP21-JP22 "A" position; JP9, JP10, JP11 "B" position.




CONSOLE COMM J8

note: only for -97 version.
 J8-5 page 4/4 E-3
 U15-20 page 3/4 B-2
 J8-9 page 4/4 E-3

SET JP12 FOR FLUOROSCAN
 +12V UNREG JP12




S	NC 12/278	F. Díaz	25/06/12	NAME	DATE	SHEET / OF	A-3024-XX								
R	NC 11/427	J.A.Garcia	25/11/11	DRAWING	F. Garcia	26/03/09	4/4	U	T	S	R	Q	P	O	← REV
Q	NC 11/036	F. Díaz	09/02/11	REVISED	A. Díaz	26/03/09									
U	NC 13/257	R. Hermosilla	19/02/14												
T	NC 13/417	R. Hermosilla	21/01/14												
REV	DESCRIPTION	ISSUED BY	DATE	SEDECAL			ATP CONSOLE								




BOARD	VERSION	JP1	JP2	JP4
3012-01	PHOTOMULTIPLIER	B	B	A
3012-02	FOR TV CAMERA	A	A	A
3012-05	4 ION CHAMBERS	C	A	B

NOTE.- JUMPER JP3 IN POSITION A FOR ION CHAMBERS WITH HIGH SENSIBILITY (GAIN 2V / mR)


AEC CONTROL SEDECAL S.A.		ENG	F. GARCIA
		CHK	A. DIAZ
DWG	REV	CN 08/137 (14/05/08)	
		A3012-01/02	REV. L
		A3012-05	REV. E
			REV.
DATE	15/05/08		



DRAWING		NAME	DATE	SHEET / OF	A3285-02	
REVISED		F. Diaz	22/02/05	1/2	E ← REV	
		A. Diaz	22/02/05			
		SEDECAL		BATTERY CHARGER		
		INNERSCAN				
REV	DESCRIPTION	ISSUED BY	DATE			



REV		DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3285-02		
					DRAWING	F. Diaz	22/02/05	2/2		
					REVISED	A. Diaz	22/02/05	E ← REV		
							BATTERY CHARGER			




90050-11


NAME		DATE		SHEET / OF		A3580-01										
DRAWING		E. Ojeda		30/05/02		1/1		Ab1	b7	b6	b5	b4	b3	b2	b1	← REV
REVISED		F. Diaz		07/07/04												
Ab1	NC 05/146	R. Frutos	20/07/05													
REV	DESCRIPTION	ISSUED BY	DATE													



BATTERY MONITOR



REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3139-01						
				DRAWING	F. Diaz	25/10/99	1/1		C	B	A	B1	← REV
				REVISED	A. Diaz	27/10/99							
C	NC 01/168	F. Diaz	21/09/01										
B	NC 00/086	F. Diaz	28/04/00										
A	NC 00/043	F. Diaz	15/03/00										
								LINE MONITOR					




90009-04

A3553-02

REV	DESCRIPTION	ISSUED BY	DATE	NAME	DATE	SHEET / OF	A3553-02				
				DRAWING	J.A.Garcia	14/02/09					
				REVISED	F.DIAZ	30/04/03	2/2				
C	NC:09/042	J.A.Garcia	23/02/09					C	B	A	← REV
B	NC:03/148	E.OJEDA	29/07/03								
A	NC:03/074	E.OJEDA	30/04/03								




SERIAL CONSOLE CONTROL



NOTE.- SEVEN SEGMENT DISPLAYS TYPE HDSP3600








		NAME	DATE	SHEET / OF	A3585-21	
		DRAWING	E.OJEDA	11/11/03	1/3	
		REVISED	A.DIAZ	11/11/03	A ← REV	
A	NC: 12/210	J.A.Garcia	25/05/12	SEDECAL		
REV	DESCRIPTION	ISSUED BY	DATE			



Conector OVERLAY

NOTE 1.- ALL LED's HLMP-0400 UNLESS SPECIFIED
 NOTE 2.- ALL DIODES 1N4148 UNLESS SPECIFIED

		NAME	DATE	SHEET / OF	A3585-21	
		DRAWING	E.OJEDA	11/11/03	2/3	
		REVISED	A.DIAZ	11/11/03	A ← REV	
A		NC: 12/210	J.A.Garcia	25/05/12	SEDECAL CONSOLE DISPLAY	
REV	DESCRIPTION	ISSUED BY	DATE			



OSCILADOR BUZZER

				NAME	DATE	SHEET / OF	A3585-21	
				DRAWING	E.OJEDA	11/11/03	3/3	
				REVISED	A.DIAZ	11/11/03		
						CONSOLE DISPLAY		
A	NC: 12/210	J.A.Garcia	25/05/12					
REV	DESCRIPTION	ISSUED BY	DATE					