




FDA Fluoroscopy Safety Standards



Hey there Subscriber,

The **FDA requires fluoroscopic systems** to include features that help operators track and control radiation exposure during procedures.

Systems must display **air kerma rates and cumulative dose** so operators can monitor patient exposure in real time. **Timers and audible signals** alert users to exposure duration, while **Last Image Hold (LIH)** and **dose display features** allow quick review without additional radiation.

These safeguards help ensure **safe, controlled, and compliant fluoroscopic imaging**.

👉 See the document for more details on FDA monitoring standards.

[View FDA Document-- 21 CFR 1020.32](#)



Have More Questions?

Browse our authorized dealer portal website. This site contains documentation to assist with installation, maintenance, and troubleshooting of Summit supported X-ray equipment.

Congratulations to the Field
Service Engineers who have
completed certifications:

Daniel Welch on Nova OTC
& SHFR

LOCATE US

SUMMIT INDUSTRIES
7555 N CALDWELL AVE
NILES, ILLINOIS 60714

CONTACT US

tscalls@summitindustries.net
customerservice@summitindustries.net
[\(773\)588-2444](tel:(773)588-2444) or [\(773\)353-4000](tel:(773)353-4000)

