




Wallstand Load Configuration

FRS WALLSTANDS NOW CONFIGURABLE ON-SITE



← Left-Hand Load



Right-Hand Load →

Hey there Subscriber,

Did you know the **FRS Wallstand** can now be configured for **either left-hand or right-hand load during installation?**

Previously, wallstands were shipped **pre-configured for a specific load direction**. Now, a single unit can be **configured on-site** to best fit the room layout.

⚠ Important: The wallstand must be **configured before installing the receptor mount**.

Please refer to **Service Manual Section 2.9.1** for instructions on configuring the wallstand for the desired load direction.



Have More Questions?

Browse our authorized dealer portal website. This site contains documentation to assist with installation, maintenance, and troubleshooting of Summit supported X-ray equipment.

SUMMITXRAY.COM →

Congratulations to the
Field Service Engineers

Scott Byers on DFMT Elite

who have completed
certifications:

LOCATE US

SUMMIT INDUSTRIES
7555 N CALDWELL AVE
NILES, ILLINOIS 60714

CONTACT US

tscalls@summitindustries.net
customerservice@summitindustries.net
[\(773\)588-2444](tel:(773)588-2444) or [\(773\)353-4000](tel:(773)353-4000)

